

Technical data sheet

Wireless accessories

Full discretion

Oticon Remote Control 3.0 offers an intuitive, straightforward and discreet way to operate hearing aids.

With the remote control, the hearing aid user can perform the most common everyday adjustments without drawing attention to the hearing aids:

- Volume up/down
- Select program
- Mute hearing aids

Full discretion is offered, since the remote control can operate the hearing aids from a distance of up to 1.8 metres without a direct line of sight between the remote control and the hearing aids.

Easy to operate

The remote control buttons are clearly visible and easy to operate.

Every time a button is pressed, a notification is heard in the hearing aids. At the same time the remote control confirms the action with a visual indication and a discreet click that can be heard and felt.

Remote Control 3.0



Remote Control 3.0 is a discreet and elegant solution for all users of Oticon hearing aids with 2.4 GHz radio technology.



Included

- Remote Control 3.0
- Batteries (2 x AAAA)
- Instructions for Use

Technical Data	
Model	RC-3A
Operating conditions	Temperature: +5 °C to +40 °C Humidity: 15% to 93 %, non-condensing Operating atmospheric pressure: 700 hPa to 1060 hPa
Storage and transportation conditions	Temperature: -25 °C to +70 °C Humidity: 15% -93%, non-condensing
Mechanical dimensions	75.3 x 35 x 12 mm; weight 32.7 g (incl. batteries)
Compatibility	Oticon hearing aids with 2.4 GHz radio technology
Battery type	2 x AAAA
Battery life at normal use	12 months
Typical transmission range	Up to 1.8 m
Status indicator	Green LED
Frequency	Licence free 2.4 GHz (ISM band)
Emission power	Max. 10 dBm EIRP



Manufactured by :
SBO Hearing A/S
Kongebakken 9
DK-2765 Smørum
Danemark
www.sbohearing.com

Imported and
Distributed by :
Audmet Canada Ltd
1600-4950 Yonge St
Toronto, ON M2N 6K1
www.oticon.ca

oticon
PEOPLE FIRST